



## High Efficiency Rust Remover

**ZAPRUST™** is a ready to use high efficiency rust remover with a fast action formula. It helps in removing and converting rust into iron and further protecting it from oxidation. The quick action of **ZAPRUST™** leaves the treated surface in a relatively passivated condition. **ZAPRUST™** restores original finish and increases use- life of instruments.



## Composition

Inorganic acid (>30% phosphoric acid) Non Ionic surfactants (< 5%),

## Contact Time

Undiluted: 5 to 10 minutes

Diluted: 10 to 30 minutes

## Features & Benefits

**ZAPRUST™** has rapid action with minimal contact time

**ZAPRUST™** is non - flammable - safe

**ZAPRUST™** doesn't contain harmful VOCs - safe to human & environment

**ZAPRUST™** is indigenously manufactured under GMP\*\* gives regular availability

**ZAPRUST™** is cost effective

## Applications

Removal of Rust from

Surgical / Medical Instruments

Dental Instruments

Industrial Equipments and Surfaces

## Usage Directions

### DILUTED

Dilute 10 % **ZAPRUST™** with potable water.

Pour **ZAPRUST™** solution in any non - metallic container.

Soak rusted part/s in **ZAPRUST™** for 10 – 30 minutes depending on surface area.

Brush to remove debris.

Thoroughly rinse with water and wipe with dry cloth or air dry.

### UNDILUTED

Pour undiluted **ZAPRUST™** solution in any non - metallic container.

Soak rusted part/s in **ZAPRUST™** for 5 – 10 minutes depending on surface area.

Brush to remove debris.

Thoroughly rinse with water and wipe with dry cloth or air dry.

**Note:** Do not soak beyond maximum recommended time.

Cat. No.	Pack size
ZAP1000	1000 mL
ZAP5000	5000 mL

