

ACCESSORIES



DISPENSER PUMPS

Disposable pumps for obtaining products from **BIOSHIELDS**® 250 mL, 500 mL & 100 mL bottles. Delivers 1.5 – 2 mL / stroke.

Cat No.: DP1001



DISPENSER STAND

Dispenser stand compatible with **STERIMAX**®, **ECOMAX**™, **PURELLIUM**™ GEL, **TRIOSEPT**™, **ALCONOX**®, **HITMAX**® and **BIOSCRUB**™ 500 mL bottles made out of plastic. Light, Robust and High Impact Resistant.

Cat No.: DS1001




IMMERSION TRAY

Immersion Tray made up of non-reactive High impact plastic. Suitable for immersion of all kinds of instruments in Hospitals, Clinics & Labs. Practical Basket Insert for easy retrieval of instruments post soak cycles. Recommended for **ENDOMAX**™, **ENDOMAX**™ 2.0, **OPAHYDE**™, **ACITAR**™ and **NOVACIDE**®

Dimensions: 72.5 cms X 23.5 cms X 10.5 cms.

Cat No.: IT1001



TRANSFER SIPHON

Reusable siphon for transferring **BIOSHIELDS**® products from refill packs to dispenser bottle packs.

Cat no.: TS1001