



## Touch - Free Liquid/Gel Dispenser

**BIOSPRAY™** is a state of art, touch-free and wall mounted dispenser to dispense handrub / handwash in medical and industrial settings.

**BIOSPRAY™** automatically dispenses both liquids and gels at a prefixed dose. This ensures adequate disinfection of hands without contaminating the environment.



## BIOSPRAY™ product description

### Features

Touch-free  
One year warranty  
After sales service  
ABS plastic  
Fixed dose dispensing  
AC adapter provided\*  
Compatible with liquids and gels

### Benefits

Prevents cross contamination  
Highly reliable  
Peace of mind  
Rust free, Durable and easily cleanable  
Adequate disinfection  
Reduced wastage of handrub / handwash  
No need of battery  
Versatile

## Specifications

Product Name : **BIOSPRAY™**  
CAT No. : BSTF 12001  
Unit Size : 142 mm (W) x143 mm (D) x255 mm (H)  
Weight : 750g (excluding chemicals, accessories etc.)  
Material : Front cover & body - ABS  
Bottle Volume : 1 Liter  
Installation : Wall mount  
AC transformer : Input : AC 100 - 240 V / Output : DC 6V

\* Battery option also available.

## BIOSPRAY™ is compatible with

### ALCONOX®

Colourless & odourless alcoholic handrub with moisturiser

### ECOMAX™

Alcoholic handrub with moisturiser

### PURELLIUM™ GEL

Alcoholic handrub gel with moisturiser

### STERIMAX®

Liquid handrub antiseptic with triple action

### TRIOSEPT™

Colourless & odourless liquid handrub antiseptic with triple action

### BIOSCRUB™

Antiseptic surgical scrub

### HITMAX®

Liquid microbicidal handwash soap

