

NOCONTA™

Air and Surface Disinfectant



Two Years' Shelf Life:

In use bottle validated for 24 months

Adjustable Spray Pattern:

Confirms optimum surface coverage

Double Sterility Assurance:

0.2 micron filtered & Gamma Irradiated

NOCONTA™

Gamma-irradiated alcoholic, trigger spray Bag in Bottle, air and surface disinfectant.

A revolving nozzle with changeable settings to provide the required spray pattern

A trigger with an ergonomic design that makes it simple to use and consistent for all operators

The outer bottle has a venting hole drilled into it so that air passes freely into the gap between the flexible bag and the outer bottle, causing the inner bag to collapse when liquid is dispensed.

Gamma Indicator

An inner flexible bag

Composition

70 % v/v Isopropyl Alcohol (2-Propanol) IP

Biocidal Activity

Broad spectrum: Bactericidal & Fungicidal

Contact Time

60 Seconds- Invisibly soiled areas
15 minutes- Visibly soiled areas

Firm outer bottle to provide structural stability.

An inner flexible bag protects its contents and creates a post-dispensing vacuum in which remaining liquid is stored.

A Division of Tulip Diagnostics (P) Ltd.



Applying Science In Disinfection

Regd. Office : Gitanjali, Tulip Block, Dr. Antonio do Rego Bagh, Alto Santacruz, Bambolim Complex P.O., Goa - 403 202, INDIA.

Website : www.bioshields.in