



# TREASURES

A brand is a voice & a product is a souvenir

\* FDA Approved  
\* Tested as per USP  
\* COA available with product

ENVIRONMENT & SURFACE  
DISINFECTANT

## AEROSEPT™

Alcohol & QAC Based Instant Surface Disinfectant  
(BKC 0.5% w/v & IPA 70% v/v)

- ▶ Safe to use on non-porous hard surfaces
- ▶ No build-up on the surface or streaking
- ▶ Contributing to overall faster contact time
- ▶ Rapid action against viruses, bacteria, & fungi



INSTRUMENT DISINFECTANT

## Acitar™

High Level Surgical Instrument Sterilizing and Disinfectant Solution  
(GA 2% w/v & BKC 5% w/v)

3 in 1 solution

- ▶ Medical instruments sterilization
- ▶ Surface disinfection
- ▶ Fogging



BODY SCRUB

## BIOSCRUB™

Antiseptic Surgical Scrub Solution  
(CHG 4% w/v & Non-Ionic Surfactants)

- ▶ CHG is recommended by CDC, USA
- ▶ Residual activity upto 6 hours
- ▶ Easy to wash off



CLEANING AGENT

## CLENZYME™

Multi-Enzymatic Cleaner for Surgical & Dental Instruments

(Enzymes 12% include Proteases, Lipases & Amylases)

- ▶ Rapidly clean surgical instruments including rigid and flexible endoscopes before sterilization
- ▶ Rapidly remove all kinds of body fluids sticking to instrument surfaces



ENVIRONMENT & SURFACE  
DISINFECTANT

## DISFIS™

Quaternary Ammonium Compound Based Surface & Environment Disinfectant  
(ADBAC 2.37% w/v & ADEAC 2.37% w/v)

- ▶ Aerial fogging in critical areas such as operation theatres (OT) and QC/QA laboratories
- ▶ Surface Disinfection
- ▶ Fogging



HANDRUB

## ECOMAX™

Alcoholic Handrub with Moisturizer  
(CHG 2.5% v/v & 70% Ethyl Alcohol)

- ▶ Broad-spectrum antimicrobial efficacy
- ▶ Synergistic power
- ▶ Long-lasting residual activity



INSTRUMENT DISINFECTANT

## Endomax™

GA Solution for High-Level Disinfection  
(Glutaraldehyde 2.45% w/v)

Low purification index ensures:

- ▶ High purity
- ▶ Efficacy
- ▶ Stability



HANDWASH



## Hitmax™

Liquid Microbicidal Handwash Soap  
(Triclosan 0.5% w/v)

- ▶ Effective cleanser for organic matter
- ▶ Optimized viscosity gives better penetration
- ▶ Powerful biocidal action



LINEN  
DISINFECTANT



Disinfectant Cleaner for Linen (DDAC 5% w/v)

- ▶ Good cleansing power
- ▶ Rapid bactericidal and fungicidal action
- ▶ Remove odor, leaving behind a pleasant smell



ENVIRONMENT & SURFACE  
DISINFECTANT

## MICROLYSE™

Disinfectant Cleaner for Floor and Hard Surfaces  
(BKC 4% w/v)

- ▶ Pleasantly perfumed aldehyde-free solution
- ▶ Effective floor disinfectant and non-staining
- ▶ Effective against bacteria and fungi



A Division of Tulip Diagnostics (P) Ltd.



Applying Science In Disinfection

INSTRUMENT DISINFECTANT



**Aldehyde-Free Cold Sterilant**  
(PHMB 3% w/v & DDAC 10% w/v)

**3 in 1 solution**

- ▶ Sterilization of surgical instruments including rigid & flexible endoscopes
- ▶ Critical area surface disinfection
- ▶ Fogging



ANTISEPTICS & SKIN PREPARATIONS



**Microbicidal Antiseptic Solution** (PHMB 1% v/v)

- ▶ Burn and Chronic wound management
- ▶ Intra-operative irrigation
- ▶ Pre & post surgery skin and mucous disinfection
- ▶ Surgical and non-surgical wound dressings
- ▶ Routine antiseptics during minor incisions catheterization, scopy etc.



INSTRUMENT DISINFECTANT



**OPA Solution for High-Level Disinfection**  
(Ortho-Phthalaldehyde 0.55% w/v)

- ▶ No activation required
- ▶ Quick turn around time
- ▶ Effective against HIV, HBV, MTB & Spores



HANDRUB



**Alcoholic Handrub with Moisturizer**

(Ethanol 62% v/v denatured with 2-propanol 3.3% v/v)

- ▶ Soft feel
- ▶ Broad-spectrum antimicrobial efficacy
- ▶ Instant biocidal action



ENVIRONMENT & SURFACE DISINFECTANT



**Aerial Fumigant, Surface, and Water Disinfectant**  
(Silver Nitrate 0.01% w/v & Hydrogen Peroxide 10% w/v)

- ▶ Aerial fogging in critical areas such as operation theatres (OT) and QC/QA laboratories
- ▶ Surface Disinfection
- ▶ Water disinfection of water reservoirs



HANDRUB



**Liquid Handrub Antiseptic with Triple Action**

(CHG 2.5% v/v, Triclosan 0.5% w/v, 2-Propanol 50% v/v, 1-Propanol 25% v/v)

- ▶ Broad spectrum antimicrobial efficacy
- ▶ No know resistance
- ▶ Synergistic power
- ▶ Long-lasting residual activity



ANTISEPTICS & SKIN PREPARATIONS



**CHX Antiseptic Solution for Preparation and Post-Closure of Skin**

(CHX 2.5% v/v & IPA 63% v/v)

- ▶ CDC, USA recommended formula
- ▶ Powerful synergistic action
- ▶ Skin safe



ENVIRONMENT & SURFACE DISINFECTANT



**Disinfectant for Hospitals, Hotels, Food and Pharma Industry**

(PHMB 3% w/v & DDAC 10% w/v)

- ▶ Critical area surface disinfection
- ▶ Fogging



RUST REMOVER



**High-Efficiency Rust Remover**

(Inorganic Acid & Non-Ionic Surfactants)

- ▶ Removal of rust from Surgical and Dental Instruments



ANTISEPTICS & SKIN PREPARATIONS



**Personal Care Antiseptic Solution**

(BKC 4% w/v)

- ▶ Green Dettol
- ▶ General antiseptic
- ▶ Non-Stinging



For further information contact :



A Division of Tulip Diagnostics (P) Ltd.  
A PerkinElmer Company

Gitanjali, Tulip Block, Dr. Antonio do Rego Bagh, Alto Santacruz,  
Bambolim Complex Post Office, Goa - 403 202, INDIA.

Tel.: +91 832 2458546-50 Fax : +91 832 2458544

E-mail : sales@tulipgroup.com Website : www.tulipgroup.com

ISO 13485  
Certified Company